

Cordless Pads



FEATURES

- Simple operation and maintenance
- For use with cura¹ Cordless Pad monitors
- No more cords between monitor and pads
- Easy to clean
- Virtually no false alarms over the lifetime of the pad
- All bed/chair pads made from non-allergenic soft PVC
- Vent holes in the neck reduce the risk of moisture ingress

TECHNICAL SPECIFICATIONS

- Cordless Chair Pad – 1 year (380 x 255)
- Cordless Bed Pad – 1 year (760 x 510)
- Cordless Pad Transmitter Insert has 50m range
- Cordless Pad Transmitter Insert fits inside cordless pads
- Cordless Pad Transmitter Insert 433MHz operating frequency
- Cordless Pad Transmitter Insert has 1 year battery life
- Cordless Pad Transmitter Insert incorporates status/in-use LED indication
- Cordless Pad Transmitter Insert incorporates on/off switch

ACCESSORIES

2718 Cordless Pad Monitor, **2556** Care Station, **2600** 30 Channel Display Panel, **2586** Hallway Alarm