

# Standard Falls Monitor with Safetime



PART NUMBER: 2701

## FEATURES

- Modern ergonomic design.
- Pause mode to attend to a patient without triggering alarm.
- Tamper-proof reset option using magnetic caregiver key.
- Choice of tone, and volume control for alarm.
- 0, 2s, and 15min. alarm delay settings.
- Port to connect with nurse call system (NC/NO).
- Operates on batteries or AC adaptor.
- Low battery notification.
- Input lost notification (Pad Lost, nurse call, AC lost).
- Change Pad (SafeTime) indication.
- Comes with mounting bracket, hanger strap, soft cover, reset key and batteries.

## TECHNICAL SPECIFICATIONS

- Power supply: 6V DC, 300mA &/or 3 x AA alkaline batteries.
- Battery life: 1 year (nominal)
- Operating temp: 5°C to 40°C
- Rel. humidity: <85%
- Housing material: ABS (UL-94HB, flammability grade)
- Dimensions (LxWxD): 120mm x 77mm x 34mm
- Net weight: 140g (without batteries), 200g (with cover)
- Pad / Mat socket: RJ10
- Adjustable volume: up to 108db @ 10cm
- Incontinence (Wetness): Under development

## ACCESSORIES

- 3516** Premium Under Mattress Pressure Pad, **2626** Premium Bed Pad, **2628** Premium Chair Pad, **2633** Premium Floor Mat, **2623** SafeTime Bed Pad - 1 Year, **2622** SafeTime Bed Pad - 90 Days, **2621** SafeTime Wheelchair / Chair Pad - 1 Year, **2620** SafeTime Wheelchair / Chair Pad - 90 Days, **2406** Plug Pack 6Vdc, **2642** Caregiver Reset Key, **3504** ActiveCare Mounting Bracket, **2636** Wheelchair / Chair Seat Belt, **2643** Nurse Call Cable.

cura<sup>1</sup>.com

cura<sup>1</sup> products are designed to augment well designed care programs and practices, cura<sup>1</sup> does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision, and are at risk.