

BRO 199 Transroller




This unique **rollator** and **wheelchair** hybrid is a clever solution for users with limitations in walking distances. The transroller can quickly and easily be converted from a rollator to a wheelchair without the need for tools.

The handles are easily height adjustable and the large 20cm wheels are suitable indoors or outdoors.

To convert to an attendant propelled wheelchair, just flip the comfortable sling backrest to the other side and fold the swing away footrests into place.



The double “X” braced frame and tough woven nylon seat means the transroller can provide temporary relief from walking for users up to 136kg.



Double “X” brace frame is the same as in wheelchairs

- Lightweight and foldable for easy transport
- Large 20cm wheels suitable for indoor or outdoor use
- The footrests are height adjustable and the legrests are swing away or can be removed
- The comfortable sling backrest conforms to the user’s back shape
- The handlebars are easily height adjustable which makes it easier to ensure user and attendant comfort
- Comes with a dual purpose removable shoulder bag

	Seat width	Seat depth	Overall height	Folded width	Product weight	Warranty	S.W.L.
BRO 199	46cm	31cm	86 - 98cm	27cm	9.3kg	1 year	136kg