

# BATH TRANSFER BENCH WITH ROTATING SEAT



Bathroom  
Support  
Aid



Patient seating from side  
angle made easier with  
rotating seat

# BATH TRANSFER BENCH C-BA422



The seat rotates 360 degrees in both directions with rotating lock at 90 degree angles for convenience.



Transfer bench shown fitted over a bath with ability to slide seat for ease of access for user.

## PRODUCT FEATURES

- ✓ ALUMINIUM FRAME
- ✓ SAFE WORKING LIMIT OF 130KG WITH SAFETY BELT
- ✓ SEAT ROTATES 360 DEGREE WITH 90 DEGREE LOCK POSITIONS
- ✓ SOFT PADDED POLYURETHANE SEATING
- ✓ SUCTION RUBBER TIPS FOR SLIP RESISTANCE
- ✓ PATIENT ACCESS FROM SIDE ANGLE IS MADE EASIER BY ROTATING SEAT

## TECHNICAL SPECIFICATIONS

**Unit Weight:** 8.8kg  
**Armrest Height:** 10cm  
**Tub Clearance:** Min 40cm - Max 50cm  
**Overall Height:** Min 93cm - Max 103cm

**Seat Width:** 43cm  
**Seat Depth:** 37cm  
**Overall Width:** 80cm  
**Overall Length:** 51cm

## CONTACT OUR SALES TEAM

**WA:** (08) 9456 1661 | **VIC:** (03) 9580 1055

**SA:** (08) 7225 8720 | **QLD:** (07) 3492 0947

**NSW:** 1300 467 737

[www.hospequip.com.au](http://www.hospequip.com.au)