

BRODA Paediatric Glider



For many children, it can be challenging and overwhelming to focus on tasks. Gliding can be an integral part of a child’s sensory needs and it provides a repetitive, soothing motion that helps to calm by allowing for gently controlled movement. Children feel more in control of their environment, aiding in self-regulation, relaxation and focus. The paediatric glider is ideal for children with Autism Spectrum Disorder, developmental delay, ADHD or anxiety. The Broda Comfort Tension Seating contours to the body for an individualised fit. The soft foam padding and supportive design also create a gentle cradling effect for additional comforting.

The gliding action ensures the base remains stationary and stable. The rear caregiver lock makes it possible to disengage the gliding action.



The Unique Comfort Tension Seating system contours to the body to create an individualised fit.

- Padded vinyl seat and back
- Smooth gliding action
- Rear transport wheels
- Caregiver lock
- Comfort Tension Seating
- Provides comfort and vestibular stimulation
- Use in classrooms, multi-sensory or quiet rooms, therapy centres, hospitals or home



	Seat Width	Seat Depth	Seat Height	Back Height	Overall Width	Overall Length	S.W.L.
Medium	36cm	31cm	31cm	51cm	50cm	70cm	91kg
Large	41cm	36cm	36cm	56cm	55cm	75cm	91kg