

TILT AND RECLINE Bariatric Shower Commode



The Revive Bariatric Tilt and Recline Shower Commode combines superior comfort with a durable stainless steel frame and components.

The padded commode seat, comfort tension seating and height adjustable headrest provide enhanced patient comfort with a weight capacity up to 363kg.

The tilt and recline feature combined with the fully elevating legrest and removable armrests provides safe positioning while allowing for thorough bathing and improved ability to wash hair and access to abdominal and peri areas.

The mid-wheel design improves manoeuvrability and prevents lateral drifting for easier transport.

- Powder coated stainless steel frame
- Adjustable back recline and posterior | anterior tilt
- Curved back to accommodate gluteal shelf
- Removable, rotational armrests
- Elevating, length adjustable legrest
- Removable padded commode seat
- 4 total-lock stainless castors
- 2 x 8" locking middle wheels with emergency brake

The infinitely adjustable 65° back recline accommodates positioning needs to allow for safe and thorough bathing and facilitates access to difficult to reach areas.

Posterior tilt improves patient posture, helps redistribute pressure and aids in falls prevention.



Seat Width	Width between armrests	Seat Depth	Seat Height	Overall Width	Tilt	Recline	S.W.L.
68cm	68 - 83cm	53 58cm	53 58cm	83 98cm	-8° to 14°	65°	363kg