

HBA 402 Shower Chair



- Frame made from ultra light aluminium tubing
- High density plastic seat with perforations
- Thicker, durable powder coated finish protects against scratches
- Height adjustable extension legs with slip resistant rubber suction tips
- Frame shape maximises space between arm rests for added comfort

HBA 403 Shower Stool



- Constructed from ultra light aluminium tubing
- Thicker powder coating finish for extra durability
- Frame shape maximises space between arm rests
- Height adjustable extension legs
- Fitted with slip resistant suction cups
- Moulded seat with perforations for drainage

	Seat width	Overall width	Seat depth	Seat height	Product Weight	Warranty	S.W.L.
HBA 402	49cm	60cm	41cm	44 - 59cm	3.3kg	1 year	140kg
HBA 403	50cm	58cm	22cm	49 - 64cm	2.8kg	1 year	140kg