





Invacare scooters cover a wide range of mobility needs!






CONTACT US

WA  (08) 9456 1661
 saleswa@hospequip.com.au

VIC  (03) 9580 1055
 vicadmin@hospequip.com.au

 hospequip.com.au

SA  (08) 7225 8720
 sareception@hospequip.com.au

QLD  (07) 3448 9243
 vicadmin@hospequip.com.au

NSW  1300 467 737
 nswadmin@hospequip.com.au

BONUS ONE-YEAR FREE OF BREAKDOWN ASSISTANCE INCLUDED*

Every new scooter purchase receives ONE YEAR FREE breakdown assistance through 24/7 Road Service Australia!

**24/7 ROAD SERVICE AUSTRALIA**

*Through our partner Invacare, HospEquip offers 1 year of free of breakdown assistance when you purchase your new scooter. To qualify customers need to have current and active policies with 27/7 Road Service Australia.

Road assistance provides assistance only when:

- Customer has an active policy
- Brake down (mechanical only) occurs outside of home





Roadside assistant provides a towing service to bring the scooter home to the customer. Repairs are not offered and cannot offer any service when the customer is located in their home.



To be eligible please ensure you request HospEquip register your warranty at the time of purchase. Once your warranty has been registered, you will be automatically signed up for Roadside assistance.

HOSPEQUIP SERVICE & REPAIRS



For our WA & VIC Customers, HospEquip offers you the option of choosing an ongoing service contract to cover preventative maintenance and therefore ensure equipment downtime is minimised. Visit our website or call us for further details.



**CONTACT US**


WA  (08) 9456 1661
 saleswa@hospequip.com.au

VIC  (03) 9580 1055
 vicadmin@hospequip.com.au

 hospequip.com.au

SA  (08) 7225 8720
 sareception@hospequip.com.au

QLD  (07) 3448 9243
 vicadmin@hospequip.com.au

NSW  1300 467 737
 nswadmin@hospequip.com.au

INVACARE Colibri Mobility Scooter

The Colibri is a simple yet stylish lightweight scooter and is ideal for those who enjoy their independent lifestyle. With its small footprint, it is extremely manoeuvrable in tight or limited spaces and can easily be dismantled and packed into the car boot. This micro lightweight scooter is not only extremely portable but powerful enough to reach speeds up to 8 km/h.



| Seat width | Seat depth | Seat height | Overall width | Overall length | Battery | Product Weight | S.W.L |
|------------|------------|-------------|---------------|----------------|-----------|----------------|-------|
| 47cm | 40cm | 37 - 42cm | 51cm | 101cm | 2 x 18 Ah | 49kg | 136kg |



INVACARE Leo Mobility Scooter



The Leo is an affordable lightweight scooter ideal for both indoor and outdoor use. It can easily be dismantled and reassembled without tools and fits neatly into the boot of most cars. The adjustable tiller and ergonomic handlebars provide greater comfort and help prevent wrist and hand fatigue.





| Seat width | Seat depth | Seat height | Overall width | Overall length | Battery | Product Weight | S.W.L |
|------------|------------|-------------|---------------|----------------|-----------|----------------|-------|
| 47cm | 41cm | 46 - 56cm | 59cm | 122cm | 2 x 36 Ah | 83kg | 136kg |



CONTACT US


WA  (08) 9456 1661
 saleswa@hospequip.com.au

VIC  (03) 9580 1055
 vicadmin@hospequip.com.au

 hospequip.com.au

SA  (08) 7225 8720
 sareception@hospequip.com.au

QLD  (07) 3448 9243
 vicadmin@hospequip.com.au

NSW  1300 467 737
 nswadmin@hospequip.com.au

INVACARE
Pegasus Metro Mobility Scooter

The Pegasus Metro is ergonomically designed and focuses on providing a safe, easy, and enjoyable ride, whether you're heading off to the local shops or a quick trip to meet up with friends. The suspension system allows you to comfortably navigate a variety of everyday surfaces and terrains and the powerful motor safely handles obstacles or uneven ground.



| Seat width | Seat depth | Seat height | Overall width | Overall length | Battery | Product Weight | S.W.L |
|------------|------------|-------------|---------------|----------------|----------------|----------------|-------|
| 51cm | 47cm | 44 -51cm | 66cm | 127cm | 2 x 12V / 50Ah | 110kg | 136kg |



INVACARE
Pegasus Pro Mobility Scooter



The ergonomically designed Pegasus PRO is ideal for longer trips into town or afternoons driving around the park with the grandchildren. It includes all the usual comfort, safety and reliability features but with extra performance. The 12" pneumatic tyres, combined with the advanced suspension system ensure you enjoy a smooth ride.





| Seat width | Seat depth | Seat height | Overall width | Overall length | Battery | Product Weight | S.W.L |
|------------|------------|-------------|---------------|----------------|----------------|----------------|-------|
| 51cm | 47cm | 44 - 51cm | 65cm | 132cm | 2 x 12V / 75Ah | 136kg | 160kg |



CONTACT US

WA  (08) 9456 1661
 saleswa@hospequip.com.au

VIC  (03) 9580 1055
 vicadmin@hospequip.com.au

 hospequip.com.au

SA  (08) 7225 8720
 sareception@hospequip.com.au

QLD  (07) 3448 9243
 vicadmin@hospequip.com.au

NSW  1300 467 737
 nswadmin@hospequip.com.au

INVACARE Comet Alpine Mobility Scooter

The Comet Alpine Scooter is the perfect solution for those who need a robust, powerful vehicle to adventure across terrains too challenging for traditional scooters. With powerful 75Ah batteries, larger 13" tyres, a high-back ergonomic seat and soft-riding full suspension, the Comet can cover a greater distance while ensuring continuous comfort.



| Seat width | Seat depth | Seat height | Overall width | Overall length | Battery | Product Weight | S.W.L |
|------------|------------|-------------|---------------|----------------|----------|----------------|-------|
| 51cm | 47cm | 44 -51cm | 68.5cm | 145cm | 2 x 75Ah | 143kg | 160kg |



INVACARE Comet Ultra Mobility Scooter



The carefully thought-out design and robust construction using the highest quality components and safety features means safety and stability are never a concern when using the Comet Ultra, even for users weighing up to 220kg. The adjustable armrests and seat columns are reinforced to offer additional support to the user when getting on and off the scooter.





| Seat width | Seat depth | Seat height | Overall width | Overall length | Battery | Product Weight | S.W.L |
|------------|------------|-------------|---------------|----------------|----------|----------------|-------|
| 66cm | 53.5cm | 44 - 51cm | 68.5cm | 145cm | 2 x 75Ah | 148kg | 220kg |



CONTACT US

WA  (08) 9456 1661
 saleswa@hospequip.com.au

VIC  (03) 9580 1055
 vicadmin@hospequip.com.au

 hospequip.com.au

SA  (08) 7225 8720
 sareception@hospequip.com.au

QLD  (07) 3448 9243
 vicadmin@hospequip.com.au

NSW  1300 467 737
 nswadmin@hospequip.com.au