

SMARTCARE PROGRAMMABLE TURNING BED



Turning Bed

The SmartCare Programmable Turning Bed is ideal for patients at risk of pressure injuries from immobility.

The SmartCare Turning Bed provides 14 degrees of rotation left & right, with both manual and automatic program cycles that can be set at 30/60/90/120 minute cycles. The cycles deliver undetectable and continuous repositioning of the patient. The SmartCare Programmable Turning Bed allows patients to be continuously repositioned without being disturbed, staff are free from time-consuming, heavy manual handling, while the providers experience reduced costs through efficient staff practices.

Contemporary best practice for pressure injury prevention is frequent, small repositioning of the patient, rather than “side to side” every 2 hours. Pressure Injuries have a very high clinical impact on patients, consuming a significant amount of the carers time and at a very high cost to providers.

SmartCare Programmable Turning Bed Pressure Care Management



Provides 14 degree rotation left and right



Automated 30/60/90/120 minute program cycles



CONTACT US

WA ☎ (08) 9456 1661
✉ saleswa@hospequip.com.au

VIC ☎ (03) 9580 1055
✉ vicadmin@hospequip.com.au

🌐 hospequip.com.au

SA ☎ (08) 7225 8720
✉ sareception@hospequip.com.au

QLD ☎ (07) 3448 9243
✉ qldadmin@hospequip.com.au

NSW ☎ 1300 467 737
✉ nswadmin@hospequip.com.au

Although treatable if found early, without specialist medical attention Pressure Injuries can be life-threatening. In addition to assisting with managing Pressure Injuries, The SmartCare Programmable Turning Bed addresses the clinical challenges of pressure management, pain management and sleep deprivation forced upon patients who are immobile for a variety of medical reasons, including Acquired Brain Injury, Spinal Injury, Burns, Intensive Care, MS, Huntingtons, Motor Neuron Disease, or Dementia.



The SmartCare Programmable Turning Bed contributes to the healing of existing Pressure Injuries and assists in preventing new injuries occurring. By using the Smartcare Programmable Turning Bed, immobile patients are automatically repositioned, left and right continuously through 30/60/90/120-minute cycles from 14 degree left tilt to 14 degrees right tilt, or manually without being touched, even during bed linen and continence pad changes or bed baths.



Recovered sleep is the biggest patient reaction to the SmartCare Programmable Turning Bed. In residential and hospital settings the patient is not disturbed for manual repositioning. In community applications, the patient, and carers both get to remain asleep while the bed manages the repositioning.

LATERAL TILT	WIDTH - OVERALL	LENGTH - OVERALL LENGTH BASE EXTD	HEAD & FOOT TILT DOWN	HEIGHT RANGE	BACKREST	KNEE BREAK	S.W.L.
14° - 0° - 14°	102.5cm	218 237.5cm	14°	47.5 - 92cm	0° - 70°	0° - 35°	200kg

CONTACT US

WA ☎ (08) 9456 1661

✉ saleswa@hospequip.com.au

VIC ☎ (03) 9580 1055

✉ vicadmin@hospequip.com.au

🌐 hospequip.com.au

SA ☎ (08) 7225 8720

✉ sareception@hospequip.com.au

QLD ☎ (07) 3448 9243

✉ qldadmin@hospequip.com.au

NSW ☎ 1300 467 737

✉ nswadmin@hospequip.com.au