

Inovys II



The principle behind the Inovys II wheelchair is that the wheelchair must be adapted to you, not you to the wheelchair. That is why the Inovys has a modular design with numerous accessories and interchangeable options.

The armrest length and height can be adjusted as well as the legrest width. The new backrest attachment system allows for depth, height and angle adjustment.

The backrest reclines and the seat can be tilted in place to optimise comfort and reduce the risk of pressure points.

The push handles, incorporating a one piece push bar, can be raised or lowered to suit the assistant and are fitted with easy to apply hand brakes and tilt levers.

Dartex back and seat for optimal comfort

Adjustable components for easy customising:

- armrest length | height
- seat depth
- legrest width
- headrest height | depth | width
- backrest depth | height | angle
- push handle height

Optional Accessories include:

- table | safety belt | spreader | padded calf support

The tilt and recline features allow positioning options which maximise comfort and reduce the risk of pressure injury



	Seat Width	Seat Depth	Seat Height	Overall Height	Product Weight	Tilt	Recline	S.W.L.
Inovys II	40 45 50cm	46 - 52cm	46 - 53cm	102 - 120cm	36kg	0° - 21°	0° - 54°	135kg