

# G - 2

*For the semi and non-ambulant user  
High posture and pressure care needs*

Underpinned by **scientific research** and ergonomic design principles, the G-2 is a groundbreaking model which achieves new levels of **comfort** and **immersion**.

This elegant chair has been designed to ensure all key adjustments are **simple** and **effective**, meeting the needs of both users and carers. Perfect for those seated for lengthy periods and who require support and positioning assistance.



### Regulated Motion

(optional with powered version)

Automated tilt repositioning between 10°- 25° ensuring dynamic movement and continuous changing of pressure hotspots



Tilt-in-Space and Legrest Elevation - manual | electrically powered

Articulating supportive backrest

Helical motion lateral supports

Arcing legrest motion - maintains footboard positioning



Adjustable seat length: 380mm - 480mm

Adjustable footboard height and angle

Includes Dyna-Tek Intelligent air cushion

Optional leg-laterals

Weight limit - 180kg



Tilt-in-space	Seat Widths	Seat Depths	Seat Height	Footboard Height	Back Height	Overall Width	Weight Capacity
50°	400   450mm 500   550mm	380 - 480mm	530mm	380 - 540mm	800mm	700mm 800mm	180kg