

You can use the following **SELECTION MATRIX** of features to help you find the right dynamic chair... though a full assessment is always recommended.

Some features are common to all the chairs in the range:

Modular design - refurbish and replace over time

Postural support and pressure relief

Simple to use | Easy to clean

Versatility - domestic use | multi-user environments | clinical settings | care settings

	Florien II	Florien Elite	G - 2	Encora	Omega	Duo
Removable side	✓	✓				
Footrest height adjustment	✓	✓	✓	✓		✓
Legrest elevation	✓	✓	✓	✓		✓
Integral pressure relief	✓	✓	✓	✓	✓	✓
Tilt-in-space - Manual	✓	✓	✓	✓		✓
Tilt-in-space - Electric			✓	✓		✓
Regulated Motion			✓			✓
Back rest recline		✓		✓		✓
Seat height adjustment		✓		✓		Electric Only
Seat depth adjustment	Depth options available	✓	✓	✓	Depth options available	✓
Seat width adjustment	Width options available	✓	Width options available	Width options available	Width options available	✓
Back height adjustment					At time of manufacture	✓
Lateral Support	Optional	Optional	Built-in	Built-in	Built-in	Optional
Waterfall back option	✓	✓		✓		✓
Weight capacity	180kg	160kg	180kg	160kg	150kg	160kg - Major 160kg - Minor 57kg - Mini

* Other fabrics available but this will result in longer delivery times and possible increased costs

To ensure the pressure relieving properties of your chair are working to their full potential, we prefer Dartex on all user **contact surfaces**.

Dartex is a strong, multi stretch, vapour permeable, waterproof fabric which is simple to wipe clean and truly complements the pressure relief incorporated into the chairs.

