


# Pads & Mats - Cordless

**cura<sup>1</sup>**  
Helping To Care



***... falls prevention for frail, elderly & dementia patients***

# Pads & Mats - Cordless

Pressure Pads & Mats make life safer for the frail, elderly & anyone at risk of falls & wandering, by alerting carers when they exit a bed, chair or step on a mat.

## CORDLESS PADS & MATS



**2629 Cordless Bed Pad**



**2630 Cordless Wheelchair / Chair Pad**




**3515 Cordless Under Mattress Pressure Pad**



**2624 Cordless Heavy-Duty Non-Slip Floor Mat**



**2719 Cordless Pad / Mat Transmitter**



*Full width bed pressure pads means reliable operation with practically no false alarms.*

## CORDLESS PAD & MAT BENEFITS

- No more cords between pad / mat & monitor! Provides a more discreet & professional appearance with less tripping hazards.

### Bed & Chair Pads

- Triggers an alarm when patient gets out of bed or chair.
- Durable & long-lasting.
- Made from non-allergenic soft PVC.
- Care instructions printed on pads.
- Vent holes in the neck reduce the risk of moisture ingress.

### Floor Mat

- Triggers alarm when patient steps on mat.
- Durable & long-lasting.
- Surface material foam PU (polyurethane).
- Bottom material PVC without phthalates.
- Wheelchair friendly, low profile edges.


## FOR MORE PADS & MATS ALSO SEE:



Universal Range



Premium Range



Economy Range



2670 PreAlert Cordless Bed Pad



SafeTime Range

## COMPATIBLE PRODUCTS

- 2700 Six Channel Cordless Falls Monitor
- 2702 Three Channel Cordless Falls Monitor
- 3501 ActiveCare LCD Pager
- 2508 LED Pager
- And all cura<sup>1</sup> cordless receivers.



Certifications may vary between products.