

Pegasus Metro Scooter



The Pegasus Metro is ergonomically designed and focuses on providing a safe, easy and enjoyable ride, whether you're heading off to the local shops or a quick trip to meet up with friends.

The suspension system allows you to comfortably navigate a variety of everyday surfaces and terrains and the powerful motor safely handles obstacles or uneven ground.

The Pegasus Metro also features automatic electronic speed reduction around bends for safer manoeuvring and the two-step disengaging lever prevents the scooter from accidentally being knocked out of drive mode and free-wheeling.

- Ergonomic steering
- Easily adjustable tiller - lever operated
- Hand brakes and clear bright brake lights
- Auto speed reduction around bends - added safety
- 11" Pneumatic tyres
- Electronics and motor protected from corrosion, dirt and water
- Two step disengaging lever prevents accidental free-wheeling
- Choice of shroud colours

Ergonomically designed steering system for easier, more responsive control and a clear easy to read LED battery level display.



Tiller height can easily be adjusted to most comfortable height to suit user.



Seat width	Seat depth	Seat height	Overall width	Overall length	Battery	Product weight	S.W.L.
51cm	47cm	44 - 51cm	66cm	127cm	2 x 12V / 50Ah	110kg	136kg