

# SOFTFORM Premier Active 2

High | Very High Risk Category



As a static mattress, the Premier Active 2 is suitable for patients considered at high risk of pressure injury development. The benefit of this system is that should it be deemed clinically necessary, the surface can be stepped up to a dynamic surface by fitting a pump to the air insert, thereby allowing additional levels of pressure relief.

This reduces a lot of moving and handling issues as the mattress can stay on the bed for a wide range of patients.

The air system is placed beneath the foam insert and the air cell inflation alternates on a 10 minute 1-in-2 cycle to provide an extremely comfortable surface while preventing pressure on any one area. The robust digital pump allows for simple plug and play activation and operates extremely quietly to limit patient disturbance. The sensory software within the pump gauges the patients' weight to maintain pressure at an appropriate level.

- Air insert within Sofiform Premier foam construction
  - step up or down between static and dynamic alternating
- New air cell design - patient comfort and reduced friction
- Mattress does not require turning
- Two way stretch PU waterproof cover
  - excellent durability and reduces shear & friction forces
- High frequency welded seams - reduces risk of fluid ingress
- Full length zip flap - conceals zip & reduces fluid ingress risk
- Toughened PU coated base - prolongs mattress longevity
- Foam "U" core - provides side wall support for patient transfer



	Width	Length	Height	Product Weight	Pump Unit Weight	Warranty	S.W.L.
Single	88cm	197cm	15cm	14kg	1.75kg	Cover - 4 years Foam - 8 years Pump - 2 years	247kg
King Single	105cm	197cm	15cm	16.7kg	1.75kg		247kg