

Titan

PATIENT WALKING TUTOR



The Titan walking tutor is the ideal solution for hospitals and institutions. The oval tube aluminium frame is extra sturdy and can withstand the constant use associated with commercial environments.

The frame is widened between the rear wheels to maximise the walking space and provide stability.

The arm supports are made of soft polyurethane with extra deep gutters for maximum comfort and safety. The press button height adjustment is simple and secure.

The hand grips can be moved forward and backwards to adjust the arm length and there are additional handles to assist the patient from a sitting position to standing.

- Heavy duty oval aluminium tube frame
- Lock brakes on 2 rear castors
- Press button directional height adjustment system
- Soft polyurethane gutter armrests
- Patient handles - assist from sitting to standing



COMBINATION ACCESSORY KIT INCLUDED:

- Storage bag
- IV pole
- Urinary drainage bag
- Vomit bag
- Oxygen tank holder



Forearm height	Length	Width	Armrest length	Castor size	Product weight	S.W.L.
95 - 128cm	84cm	70.5cm	38 - 51cm	10cm	8kg	200kg