

Utility Chair



- Frame made from sturdy steel tubing
- Powder coated for extra durability
- Height adjustable extension legs
- Seat depth and forward / rearward adjustment
- Comfortable padded seat, back and armrests
- Hard-wearing plush upholstery
- Ideal for domestic and facility use

High Back Chair



- Sturdy steel frame
- Contoured back-rest with lumbar support
- Depth adjustable seat
- Padded upholstery with wide padded armrests
- Height adjustable extending legs
- Safe working limit up to 205kg

	Overall height	Overall length	Overall width	Seat width	Seat depth	Seat height	Product Weight	Armrest height	S.W.L.
Utility (45)	89 - 105cm	49.5cm	55cm	45cm	43 - 46cm	48 - 65cm	10.2kg	15cm	205kg
Utility (52)	89 - 105cm	49.5cm	61cm	52cm	43 - 46cm	48 - 65cm	11.6kg	15cm	250kg
High Back	114 - 129cm	65cm	64.5cm	49cm	50 - 55cm	47.5 - 63.5cm	13.8kg	15cm	205kg